

Items to stock a "Cool-it" Room(H4200)

Find out where patient is coming from.

If coming from CV lab you won't need cooling device pack, temperature probes, central line kit, etc.

Item	Where to Find it
• Quad Pole with 4 pumps	• Clean utility room (or dirty)
• Cooling device machine and attached cables	• Clean utility room
• (1) Large Cooling device Pack (do not open)	• Top of stock carts (clean U)
• Lifepak 20	• Clean utility room (shelf)
• Swan set up	• Clean utility room
• Nerve Stimulator (small patches in box)	• Red bins (clean U)
• (2) Narcotic Lock Boxes (may need to order)	• Upper red bins (clean U)
• Doppler and gel	• Lower red bins (clean U)
• (3-4) suction regulators, canisters, wall rings	• Red bins for regulators • By wire rack for canisters
• Graduated cylinders (3) • (1) for NS syringes • (1) for 60mL syringe • (1) for safe set cannulas	• Stock cart "Personal Care"
• (1) Liter of normal saline	• In stock cart marked "IV"
• Primary & Secondary IV tubing	• In stock cart marked "IV"
• (10-12) pre-filled 10mL NS syringes	• In stock cart marked "IV"
• Nasogastric tube	• In stock cart "GI/GU"
• (60mL) syringe	• In stock cart "GI/GU"
• Esophageal/rectal temp. probe (back up)	• In stock cart "GI/GU"
• (2) ABG kits	• In stock cart "Respiratory"
• Safe set kit with safe set blunt cannulas	• Stock cart "Central Line"
• Central line box	• Top of stock carts (clean U)
• Soft restraints	• Wire rack, shelf
• Glucometer and testing supplies	• Found in alcoves/tech area
• Rainbow of tubes/transfer device	• Tech station